

# Environment of Care product portfolio



**Sani-Prime**

DISINFECTS IN 1 MINUTE

DISINFECTANT QUAT/ALCOHOL FORMULATION

KILLS 56 MICROORGANISMS

A 1 MINUTE OVERALL CONTACT TIME, INTERMEDIATE LEVEL DISINFECTANT FOR ALL YOUR HIGH-TOUCH, HOUSEWIDE SURFACE DISINFECTION NEEDS

**Sani-HyPerCide**

DISINFECTS IN 1 MINUTE

EFFECTIVE AGAINST C. DIFF 5 MINUTES

KILLS 24 MICROORGANISMS

STANDARDIZE YOUR DISINFECTION PROTOCOLS WITH THIS POWERFUL H<sub>2</sub>O<sub>2</sub> FORMULA THAT KILLS *CLOSTRIDIODES DIFFICILE* WITHOUT COMPROMISING COMPATIBILITY

**Super Sani-Cloth**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN 2 MINUTES

#1 DISINFECTANT WIPE IN HEALTHCARE

DISINFECTANT QUAT/ALCOHOL FORMULATION

KILLS 32 MICROORGANISMS

IDEAL FOR DAILY USE IN FAST-PACED ENVIRONMENTS THAT REQUIRE SHORT CONTACT TIMES AND BROAD COVERAGE OF MICROORGANISMS

**Sani-Cloth AF3**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN 3 MINUTES

DISINFECTANT QUAT ALCOHOL FREE FORMULATION

KILLS 45 MICROORGANISMS

PROTECTING YOUR PATIENTS HAS NEVER BEEN EASIER. FRAGRANCE FREE - IDEAL FOR USE AROUND PATIENTS AND STAFF, ESPECIALLY THOSE WITH RESPIRATORY SENSITIVITIES

**Sani-Cloth Plus**  
GERMICIDAL DISPOSABLE CLOTH

DISINFECTS IN 3 MINUTES

DISINFECTANT QUAT FORMULATION

KILLS 18 MICROORGANISMS

IDEAL FOR USE IN ALTERNATE CARE SETTINGS INCLUDING PHYSICIAN AND DENTAL OFFICES

**Sani-Cloth Bleach**  
GERMICIDAL DISPOSABLE WIPE

DISINFECTS IN 4 MINUTES

DISINFECTANT 1:10 BLEACH DILUTION FORMULATION

KILLS 50 MICROORGANISMS

IDEAL FOR DISINFECTING AREAS CONTAMINATED WITH *CLOSTRIDIODES DIFFICILE* SPORES AND NOROVIRUS

**Easy Screen**  
CLEANING WIPE

IDEAL FOR ANY EQUIPMENT REQUIRING 70% IPA

CLEANING 70% ISOPROPYL ALCOHOL FORMULATION

A COMPATIBLE 70% IPA TOUCHSCREEN CLEANER THAT WORKS AT THE SPEED OF TECHNOLOGY<sup>3</sup>

\* Effective against *Candida auris* (*C. auris*).

<sup>1</sup>Kills SARS-CoV-2 (Coronavirus), the virus that causes COVID-19, on hard, nonporous surfaces excluding Easy Screen® wipes; <sup>2</sup>Market leadership claims based on Q1 2021 Clarivate Data; <sup>3</sup>Sani-Cloth® Plus Germicidal Disposable Wipe compares to competitor Quat/Low-Alcohol disinfectant products; <sup>4</sup>Cleaning the screen of your mobile device with a wipe containing alcohol may reduce the lifespan of any anti-fingerprint coating. Please check your device manufacturer's cleaning instructions before using Easy Screen® wipes and all other PDI disinfectant products.

## Sani-Prime

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P25372	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P24284	7.5" x 15"	70/can 6 cans/case
Large Individual Packets	H06182	5" x 8"	50/box 10 boxes/case
X-Large Individual Packets	U13195	11.5" x 11.75"	50/box 3 boxes/case

## Sani-HyPerCide®

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P27372	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P26584	7.5" x 15"	65/can 6 cans/case

## Super Sani-Cloth®

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	Q55172	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P86984	7.5" x 15"	75/can 6 cans/case
Large Individual Packets	H04082	5" x 8"	50/box 10 boxes/case
X-Large Individual Packets	U87295	11.5" x 11.75"	50/box 3 boxes/case
Softpack	A22480	8.2" x 9.8"	80/pack 9 packs/case

## Sani-Cloth® AF3

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	P13872	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	P72584	7.5" x 15"	75/can 6 cans/case
Large Individual Packets	H59200	5" x 8"	50/box 10 boxes/case
X-Large Individual Packets	U27500	11.5" x 11.75"	50/box 3 boxes/case
Softpack	M8063S80	8.2" x 9.8"	80/pack 9 packs/case
Pail	P1450P	7.5" x 15"	160/pail 2 pails/case
Refill (for Pail)	P2450P	7.5" x 15"	160/refill 2 refills/case

## Sani-Cloth® Plus

GERMICIDAL DISPOSABLE CLOTH

	REORDER NO.	WIPE SIZE	CASE PACK
Large Canister	Q89072	6" x 6.75"	160/can 12 cans/case
Extra-Large Canister	Q85084	7.5" x 15"	65/can 6 cans/case

## Sani-Cloth® Bleach

GERMICIDAL DISPOSABLE WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Clinical Size Wipes Canister	P84172	6" x 5"	160/can 12 cans/case
Large Canister	P54072	6" x 10.5"	75/can 12 cans/case
Extra-Large Canister	P25784	7.5" x 15"	65/can 6 cans/case
Large Individual Packets	H58195	5" x 7"	40/box 10 boxes/case
X-Large Individual Packets	U26595	11.5" x 11.75"	40/box 3 boxes/case
Pail	P7007P	7.5" x 15"	160/pail 2 pails/case
Refill (for Pail)	P700RF	7.5" x 15"	160/refill 2 refills/case

## Easy Screen®

CLEANING WIPE

	REORDER NO.	WIPE SIZE	CASE PACK
Medium Canister	P03672	6" x 9"	70/can 12 cans/case

### Educational and point-of-care accessories.

Developed to assist you in facility protocol compliance.



**PDI Sani-Bracket®**  
for wall/pole mounting of a canister



**Sani-Canister Caddy®**  
for placement on flat surfaces (bilingual)



**Sani-TAG® Equipment ID System**  
for identification of equipment to ensure the proper wipe is used on a particular piece of equipment



**Educational Wall Charts and signage (bilingual)**

For more information contact your local PDI representative.

400 Chestnut Ridge Road, Woodcliff Lake, New Jersey 07677 T: 800.999.6423 W: pdihc.com  
©2023 PDI UPDATE 0423 PDI09190835